

Water Bottle System Clean Water

Adapter Kit for ADEC® Bottle System



Capacity of water bottle:
1 litre

Mytronic bottle system



ADEC® Bottle System



Easy mounting to an existing ADEC® bottle system

- **Water bottle receptacle for ADEC® bottle system**
 - easy conversion of the ADEC® bottle system to the Mytronic bottle system
 - water bottle is connected to compressed air (regulated by an automatic pressure regulator)
 - air pressure stops automatically when water bottle is exchanged

Individually configurable nipples for water and air



Mytronic water bottle receptacle for ADEC® bottle system

All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH
 Riedelbach | Technologiepark
 In den Enstewiesen 7
 61276 Weilrod | Germany
 Tel +49 (0) 6083 95989-0
 Fax+49 (0) 6083 95989-11
 info@mytronic-dental.com