

3-Way Syringe

MC 3FP/3FP-S Syringe Angled or Straight Form



! Advantages

- Very light! Weighs only 35g (without hose)
- Syringe head can be removed for easy cleaning/disinfection
- With soft grip handle to avoid slipping during use
- Quick-disconnect feature
- To be used with disposable syringe tips, Stropko Irrigator, autoclavable metal syringe tips
- CE approved

(For further details please see reverse!)

Soft Grips available in the following colours:



All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH
 Riedelbach | Technologiepark
 In den Ensterwiesen 7
 61276 Weilrod | Germany
 Tel +49 (0) 6083 95989-0
 Fax+49 (0) 6083 95989-11
 info@mytronic-dental.com

3-Way Syringe

MC 3FP/3FP-S Syringe Angled or Straight Form

Details:

- Weight **only 35 g**
 - About 2/3 less than regular metal syringes!
- Quick disconnect feature for mounting disposable syringe tips,
- Stropko Irrigator, autoclavable metal syringe tips
- Clear separation of air and water, according to
- European standard (air on the left; water on the right) or US standard (air on the right; water on the left)
- Push buttons for proper dosing of air and water
- Almost no release of heat due to use of HTC polyacetal resin
- Soft grip to avoid slipping and ensure safe use
- Service-friendly control valves
- No metal abrasion on instrument holder
- Company logo may be imprinted (upon request)
- **Optional:** Syringe head removable for easy cleaning; Quick shut-off for water and air
- CE approved



! Sleeves available for different types of instrument holders (please inform us which syringe has been used previously)



All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH
 Riedelbach | Technologiepark
 In den Ensterwiesen 7
 61276 Weilrod | Germany
 Tel +49 (0) 6083 95989-0
 Fax +49 (0) 6083 95989-11
 info@mytronic-dental.com