

# Water Bottle System

## Water Bottle System Clean Water

for mounting ON unit (via angle) or IN unit (mounting in unit from below)

for mounting ON type with mounting bracket stainless steel or white  
 with thread for S-Stream bottles      with standard thread for e. g. DCI® bottles



**!** For hygienic safety you can rely on Mytronic's Clean Water System

- Can be adapted to any treatment unit
- Use of clean water ensures technical safety

### All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH  
 Riedelbach | Technologiepark  
 In den Ensterwiesen 7  
 61276 Weilrod | Germany  
 Tel +49 (0) 6083 95989-0  
 Fax+49 (0) 6083 95989-11  
 info@mytronic-dental.com

# Water Bottle System

## Water Bottle System Clean Water

for mounting ON unit (via angle) or IN unit (mounting in unit from below)

### For Hygienic Safety You Can Rely on Mytronic's Clean Water System

- The system can be adapted to any treatment unit
  - Mounting on treatment unit by using a mounting bracket
  - Water bottle connected to compressed air (regulated by an automatic pressure regulator)
  - Automatic pressure release when exchanging the water bottle
  
- The use of clean water ensures technical safety
  - Mytronic's water bearing hoses are approved for drinking water use and are biocompatible.
  - Our hoses are tested in compliance with DVGW worksheet no. W 270\*. (Reduction of micro-biological growth such as bio-slime, fungi and bacteria in drinking water applications).


\* DVGW = German Association of Gas and Water

- E Stainless steel version for mounting ON type with or without mounting bracket

#### Mounting IN type

with thread for S-Stream bottles

with standard thread for e.g. DCI® bottles



**!** The Mytronic Water Bottle Systems are also available with a higher water flow rate for cup filling! Please contact us for further information!

### All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

DCI® = Registered Trademark of DCI®

Mytronic GmbH  
 Riedelbach | Technologiepark  
 In den Ensterwiesen 7  
 61276 Weilrod | Germany  
 Tel +49 (0) 6083 95989-0  
 Fax +49 (0) 6083 95989-11  
 info@mytronic-dental.com