


MedtroSmart Dental Unit

MedtroSmart Dental Unit for Air Driven Instruments

- Small robust and low-maintenance dental unit
- Includes:
 - 1 x MC 3 FP-S Syringe
 - 2 x 4 hole turbine connection (for turbines, air motors or air scaler)
 - instrument hoses including oil collector and unit side connector
 - separate water bottle system „Clean Water“
 - removable pneumatic foot switch with unit side connector
- Anodised aluminium housing
- With wall fastening or mounting bracket
- Removable instrument tray for easier cleaning
- For use only air pressure necessary
- Certified medical device
- Company logo may be laser engraved (upon request)

(Instruments not included in delivery)



All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH

Riedelbach | Technologiepark
In den Ensterwiesen 7
61276 Weilrod | Germany


Tel +49 (0) 6083 95989-0
Fax+49 (0) 6083 95989-11
info@mytronic-dental.com

MedtroSmart Dental Unit

MedtroSmart Plus Dental Unit for Air Driven Instruments

- Small robust and low-maintenance dental unit
- Includes:
 - 1 x MC 3 FP-S Syringe
 - 2 x 4 hole turbine connection (for turbines, air motors or air scaler)
 - instrument hoses including oil collector and unit side connector
 - separate water bottle system „Clean Water“
 - removable pneumatic foot switch with unit side connector
- Corian housing
- With wall fastening or mounting bracket
- Removable instrument tray for easier cleaning
- For use only air pressure necessary
- Certified medical device

(Instruments not included in delivery)



All from a single source!

Mytronic has its own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH

Riedelbach | Technologiepark
In den Ensterwiesen 7
61276 Weilrod | Germany

Tel +49 (0) 6083 95989-0

Fax+49 (0) 6083 95989-11

info@mytronic-dental.com