



MYTRONIC BOTTLE SYSTEM

for mounting on unit with / without mounting bracket
for built in unit (mounting from below)

MOUNTING ON UNIT VERSION WITH MOUNTING BRACKET STAINLESS STEEL OR WHITE



! Hygienic safety thanks to the Mytronic Bottle System

- Can be converted for any treatment unit
- Technical safety through clean water



Mytronic Bottle System

HYGIENIC SAFETY THANKS TO THE MYTRONIC BOTTLE SYSTEM

Can be converted for any dental unit:

- Mounting on the dental unit with and without a mounting bracket
- Built in the dental unit
- Connection of the water bottle to the compressed air (reduced by an automatic pressure regulator)

Technical safety

- The Mytronic Bottle System is biocompatible.
- The water carrying hoses are tested according to DIN EN 15421*.
- Approved as medical device class 1 according to MDR.
- Automatic stop of compressed air for bottle change.
- Validated reprocessing.
- Easy filling and cleaning due to large bottle neck.

* Reduction of microbiological growth such as bioslime, fungi and bacteria in the drinking water area



Labeling on customer request possible!

the bottle can be labeled by printing and laser.

BUILT IN UNIT VERSION

For built in the dental unit



! The Mytronic Bottle System is also available with a higher water flow for cup filling! For more information, please contact us!

Everything from one source!

Mytronic has its own development department, its own tool and mold making and produces with an in-house production depth of over 90%. This means that product innovations can be and a high quality standard can be guaranteed. Mytronic is certified according to DIN EN ISO 13485.

Mytronic GmbH

Riedelbach | Technologiepark
In den Ensterwiesen 7
61276 Weilrod | Germany

Tel +49 (0) 6083 95989-0
Fax +49 (0) 6083 95989-11
info@mytronic-dental.com